Vscan Air™

Indications Reference Guide

This information is meant to be reference for examples of anatomies and examinations that can be evaluated by this product. The list may not be all inclusive.

Deep Scanning Transducer (Curved Array)

Complete both shallow and deep exams with a simple flip of the 2-sided probe.

The curved array transducer on Vscan Air CL supports Black/ white (B-mode), Color (Color doppler) and Harmonic imaging modes. Vscan Air is indicated for ultrasound imaging, measurement, and analysis of the human body in clinical applications that include the following.

CLINICAL APPLICATION	ANATOMY	EVALUATION
Abdominal ultrasound (Adult/ Pediatrics)	 Gall bladder, biliary tree, common bile duct Liver Pancreas Spleen Bowel including Appendix, small bowel loops Abdominal aorta Kidneys 	 Gall stones Gall bladder inflammation (wall thickening, surrounding fluid) Biliary obstruction (duct dilatation) Hepatomegaly Fatty liver Splenomegaly Intestinal obstruction Appendicitis Peritoneal fluid Mass/cyst / Abscess Abdominal aortic aneurysm Kidney stones
Urology (Adult/ Pediatrics)	 Kidneys Ureter Urinary Bladder Uretero-vesicular junction Prostate 	 Kidney, ureteral, bladder stones Kidney length Hydronephrosis Bladder dysfunction Pre-post Bladder volume Bladder inflammation (wall and mucosal changes, calcifications) Prostate size and volume Mass/cyst Ureteral jets with color
OB-Gyn	 Uterus and endometrium Ovaries Cervix Pouch of Douglas (POD) Gestational Sac (GS) Placenta Amniotic fluid Fetus(es) 	GS location (Intra-uterine/ extra-uterine) Fetal viability/ heart motion Placenta position (including low-lying and previa) Fetal position and presentation Amniotic fluid assessment Cervical length measurement/ cervical insufficiency Fetal well-being assessment: Biophysical profile (breathing, movements, tone, amniotic fluid) Confirmation of fetal death Intrauterine device position Endometrial thickness measurement Uterine/ adnexal mass/ cyst (fibroids, cysts) Free fluid in Pouch of Douglas



Deep Scanning Transducer (Curved Array) Cont.

CLINICAL APPLICATION	ANATOMY	EVALUATION
Lung/Thoracic (Adult/ Pediatrics)	 A-lines, B-lines, E-lines Pleura Lung tissue Lung sliding Lung point	 Pneumothorax and hemothorax Pleural Effusion Lung consolidation Pneumonia/ pneumonitis Pulmonary fibrosis Pulmonary interstitial and inflammatory disorders (Ex. ILD, COPD) Acute respiratory distress syndrome
Cardiac and hemodynamic assessment (Adult/ Pediatrics*)	 Heart (atria, ventricles, valves) including pericardium – Subcostal view Inter-atrial and interventricular septum Pulmonary arteries/ veins IVC 	 Pericardial fluid LV and RV size and function Valvular regurgitations/ stenosis Volume status and responsiveness IVC size Respiratory variation
Musculoskeletal (Conventional) (Adult/ Pediatrics)	 Hip/knee/ Shoulder joints Femur Humerus/elbow Tibia/fibula Radius/ulna Muscles Ligaments Tendons Nerves 	 Fluid Cyst/mass Long bone fractures Ligament and joint integrity Tendon injuries (tendonitis, rupture/tear) Muscle tears Peripheral nerve blocks
Procedure guidance (Adult/ Pediatrics)	 Heart Lung Uterus Abdomen Thorax Bladder Nerve plexus Hip/knee /Shoulder joints 	 Fluid detection: Pericardial, Pleural, Peritoneal, Amniotic, Joints Procedures: Thoracentesis, Paracentesis, Pericardiocentesis, Amniocentesis, Arthrocentesis Foreign body visualization/ localizations Bladder catheterization Nerve blocks Biopsy Placement and monitor position of tubes and catheters
Protocols	 Heart IVC Lungs Abdomen	• FAST • eFAST • BLUE • FASH • FASE



Shallow Scanning Transducer (Linear Array)

Complete both shallow and deep exams with a simple flip of the 2-sided probe.

The linear array transducer on Vscan Air CL supports Black/ white (B-mode), Color (Color doppler) and Harmonic imaging modes. Vscan Air is indicated for ultrasound imaging, measurement, and analysis of the human body in clinical applications that include the following.

CLINICAL APPLICATION	ANATOMY	EVALUATION
Peripheral Vascular (Adult/ Pediatrics)	 Arteries including Carotid, vertebral, subclavian, axillary, brachial, iliac, saphenous, popliteal, femoral Veins including Jugular, subclavian, cephalic, basilic, saphenous, femoral, popliteal, tibial 	 Deep vein thrombosis Atherosclerosis- Intima media thickness, plaques, vessel occlusion/ stenosis Subclavian Steel syndromes
Lung/ Thoracic (Adult/ Pediatrics)	A, B, E linesPleuraLung tissueLung slidingLung point	Pneumothorax and hemothorax Pleural Effusion Lung consolidation Pneumonia/ pneumonitis Pulmonary fibrosis Pulmonary interstitial and inflammatory disorders (Ex. ILD, COPD) Acute respiratory distress syndrome
Small organs (Adult/ Pediatrics)	Testes Scrotum Thyroid Breast Bowel Abdominal wall Skin Subcutaneous tissue Fascia Lymph nodes	Testicular torsion (size, echo-texture and vascularity) Epididymo-orchitis Fluid collection in scrotal sac Hematomas, hernias Breast nodules, mass, cyst Abdominal wall masses, hernias Thyroid nodules/cyst/mass/ diffuse enlargement Bowel pathology (ex. appendicitis, diverticulitis, intestinal obstruction) Pyloric stenosis/ Intussusception for pediatric patients) Soft tissue infection (cellulitis, abscess, bed sore) Foreign body visualization/ localization) Cutaneous mass
Musculoskeletal – (Superficial and conventional) (Adult/ Pediatrics)	 Tendons Muscles Ligaments Nerves Long bones (ex. Humerus, Radious, Ulna, Femur, Tibia, Fibula) Joints (ex. Ankel, Shoulder, Knee, Elbow, Wrist) Joint space/bursae 	 Tendon injuries (tendonitis, rupture/ tear) Muscle tears Long bone fractures Carpal Tunnel syndrome Fluid collection in joint space, muscles, bursae Joint and ligaments integrity Cyst/ mass Hip joint evaluation for neonates and infants
Nerves (Adult/ Pediatrics)	Peripheral nerves including examples as Interscalene, supraclavicular, infraclavicular, axillary plexus, Median N, Radial N, Ulnar, Femoral, Popliteal, Tibial, Peroneal, Saphenous N	Peripheral nerve blocks



Shallow Scanning Transducer (Linear Array) Cont.

CLINICAL APPLICATION	ANATOMY	EVALUATION
Neck and airway (Adult/ Pediatrics)	 Cervical Lymph nodes Trachea Epiglottis, cricoid cartilage, cricothyroid membrane Esophagus Vocal folds 	Neck massessAirway assessmentVocal cord dysfunction
Procedural guidance (Adult/ Pediatrics)	 Thorax Veins (including Jugular/ Subclavian/ Axillary/ Femoral / Brachial/ Basilic/ Cephalic) Arteries (including femoral, radial, brachial, axillary, dorsalis pedis) Peripheral nerves Joints Vertebral spaces Skin and subcutaneous tissue Trachea and surrounding structures 	Fluid detection and removal support: thoracentesis Peripheral venous access Central venous catheterization Arterial access Assessment and support of dialysis access Nerve blocks Joint aspiration and injections Cyst aspiration Biopsy Abscess drainage Foreign body visualization/ localization Lumbar Puncture Endotracheal tubes placement and confirmation Support placement and monitor position of tubes and catheters
Ophthalmic**	Optic nerve sheath Retina Globe Lens	 Retinal detachment Vitreous hemorrhage Intra-ocular foreign body visualization Globe rupture Optic Nerve sheath diameter Lens dislocation
Cephalic (Neonatal)	Fontanelle Superficial and mid-superficial cranial structures	Gyral-sulcal anatomy Superior sagittal sinus thrombosis Cerebral edema Extra-axial fluid collections
Protocols	• Lungs	• eFAST • BLUE

 ${}^{\mbox{\tiny †}}\mbox{Pediatric population for Cardiac application defined}$ as minimum body weight 40 Kg and above.

www.gehealthcare.com

© 2021 General Electric Company.

All rights reserved. GE, the GE Monogram and Vscan Air are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc. doing business as GE Healthcare.

Products mentioned in the material may be subject to government regulation and may not be available in all countries. Shipment and the effective sale can only occur after approval from the regulator. Please check with your local GE Healthcare representative for details.



 $[\]ensuremath{^{^{**}}}\xspace$ Ophthalmic not available in Japan and China